



## GROWTH MEDIA

### 4 - Tryptic soy broth (TSB)

Components	Concentration (g/L)
Pancreatic digest of casein	17
Papaic digest of soybean	3
NaCl	5
K <sub>2</sub> HPO <sub>4</sub>	2.5
Glucose	2.5

#### PROCEDURE:

Dissolve 30 g of Tryptic soy broth powder or the quantities of each component listed in the table in distilled H<sub>2</sub>O, heat the mixture while stirring to fully dissolve all components and adjust the volume to 1L. Autoclave at 121 °C for 20 min.